

In solution protein digestion using trypsin as protease

Urea is a chaotropic agent and disrupts three dimensional structure of proteins and denatures them. **However, urea + heat + protein = carbamylation;** urea in solution is in equilibrium with ammonium cyanate, that may decompose to ammonia and isocyanic acid (HNCO). Isocyanic acid attach the N-terminal of the protein, but also the side chains of lysine and arginine residues rendering a protein unsuitable for many enzymatic digests ($\text{HN}=\text{C}=\text{O} + \text{H}_2\text{N}\sim \rightarrow \text{H}_2\text{N}-\text{CO}-\text{NH}\sim$). Urea will always degrade to isocyanic acid, so urea solutions must be made fresh, and it is recommended to add 20mM methylamine (CH_3NH_2) to the urea solution prior to use (urea can also be removed before digestion using reversed phase chromatography)

Protein solvation/denaturation (applies for 100 µg protein or lower)

Dissolving the protein pellet;

The pellet may be difficult to dissolve. Add **20µl urea solution** (see right panel) and pipette gently up and down, sonicate if necessary.

Add **20µl trypsin buffer** (see right panel), and incubate at RT in Eppendorf mixer for 5 min (slow agitation).

Urea solution; 8M Urea/20mM methylamine;

Add **480 mg Urea** (art. no. 51458, Sigma-Aldrich), **1.7µl 40 wt% methylamine in H₂O** (art. no. 426466, Sigma-Aldrich) and **630µl dH₂O**.

Trypsin buffer; 50mM Tris/1mM CaCl₂;

Add **0.61g Tris** (art. no. 252859, Sigma-Aldrich) and **15mg CaCl₂ x 2H₂O** (art. no. 21097, Sigma-Aldrich, inhibits chymotrypsin activity) to about 90ml dH₂O. Correct the pH to 7.8-8 with HCl and adjust the volume to 100ml. Store the solution at 4 °C.

Reduction and alkylation

Reduction;

Add **4µl 100 mM DTT** (*see right panel*), and incubate for 1 hour at room temperature (do NOT use 56°C as with gel pieces. That will cause carbamylation due to the presence of urea in the sample).

100 mM DTT in MilliQ water:

Add **15.4 mg DTT** (*DiThioThreitol, art. no. 171318-02, Amersham Biosciences*) to 1ml dH₂O (*may be aliquoted as a 1M solution, and kept in freezer*).

Alkylation;

Add **5µl 200 mM IAA** (*see right panel*) for cystein alkylation, and incubate for 1 h at room temperature (dark).

200 mM IAA in MilliQ water:

Add **18.5mg IAA** (*Iodoacetamide, art. no. I-6125, Sigma Aldrich*) to 0.5ml dH₂O (*must be freshly made and kept in the dark*).

To avoid unwanted protease alkylation, add 0.8µl 100 mM DTT, and incubate 10 min. at room temperature.

Digestion

Sample dilution;

Add 110.2µl Trypsin buffer (the urea concentration is now 1M).

Trypsin;

Add trypsin at a concentration about 50 times lower than the amount of protein in the sample. If the sample contains approx. 100 µg protein, add 2µg of protease (*see right panel*). Measure pH using an indicator paper (litmus paper or similar), and incubate samples at 37°C overnight on a shaker

2µg Trypsin Porcine (4µl) (Promega, art. no. V 5111);

Dissolve each ampoule (20 µg trypsin porcine) in 40 µl 50 mM acetic acid (*resuspension buffer supplied from Promega with the trypsin powder*). The trypsin concentration in this stock solution is then 0.5 µg/µl

Acidification

In this final step, add 15 µl 10% FA (formic acid) to quench the digestion activity. We now have approximately 0.5 mg/ml digested protein solution at pH 3. The Urea concentration in this solution (below 1M) allows analysis directly by MALDI or LC-MS. The solution should be desalted/concentrated on reversed phase microcolumns before either MALDI-ToF or nanoflow LC-MS.